



FRESENIUS
KABI



freeflex[®]+


with Luer Lock Port

Together, We Can Care Better



freeflex+ Swabbable Luer Lock Port

- + Luer Lock-compatible injection port helps prevent needle stick injuries and minimizes risk of bag damage, leakages and contamination of patients or staff
- + Swabbable, reusable injection port helps ensure safe and convenient disinfection between multiple drug additions
- + Self-sealing septums in both ports help prevent leakages
- + Port accommodates a variety of male Luer Lock syringes



Easy Push & Turn Attachment

- + Only half a turn of a Luer Lock-equipped transfer adapter or syringe is needed to reach the end position for a safe and tight connection with the bag
- + Minimal force is required to push through the self-sealing septum with a syringe or vial adapter



freeflex®+ IV bags are designed to help reduce the risk of infusion errors, injury and waste. **So together, we can care better.**

- + Needleless connection** – Luer Lock port helps reduce the risk of needle stick injury or bag puncture.
- + Self-sealing septums** – Both **freeflex+** IV ports contain a self-sealing septum to help prevent leaks after removing the Luer Lock or spike, potentially reducing risk of hazardous exposure.
- + Swabbable, reusable Luer Lock port** – Easy to disinfect between injections. Up to 5 injections possible without affecting the resealing property of the injection port.
- + Sterility protection** – Both ports contain easy break-off caps, providing visual tamper evidence to help ensure sterility prior to use.



Helping improve the safety of IV drug delivery.

freeflex+ IV Bags are not made with PVC or DEHP¹

- Minimizes patient exposure to toxic DEHP, contained in PVC¹



freeflex+ IV bags can be used across your facility for a broad clinical application.

1. Healthcare Without Harm. Why Health Care is Moving Away from Hazardous Plastic Polyvinyl Chloride (PVC). 6 April 2006.

We asked pharmacists and nurses what they wanted in an IV bag. **Then we built it.**

To design **freeflex+** IV bags, Fresenius Kabi collaborated with nurses, pharmacists and pharmacy technicians. We listened carefully to their issues and frustrations regarding their day-to-day experiences. Together, we designed simplicity and safety into every detail of the innovative **freeflex+** IV bag.

All the features and advantages of freeflex – with a unique Luer Lock injection port to facilitate needleless drug addition. A plus in safety for you and your patients.



Added Features of freeflex+

- Unique light blue and dark blue ports indicate flow and function to help prevent errors and have wide spacing for easy handling.
- Clear, easy-open overwrap.

freeflex+ Product Portfolio

0.9% Sodium Chloride Injection, USP



Product Code	NDC Code	Fill Size (mL)	Units per Case
232150	65219-218-50	50	60
232200	65219-220-10	100	50
232225	65219-222-25	250	30
232250	65219-224-50	500	20

0.45% Sodium Chloride Injection, USP



Product Code	NDC Code	Fill Size (mL)	Units per Case
242250	65219-226-50	50	60
242200	65219-228-10	100	50
242325	65219-230-25	250	30
242350	65219-232-50	500	20

5% Dextrose Injection, USP



Product Code	NDC Code	Fill Size (mL)	Units per Case
252300	65219-234-10	50	60
262300	65219-236-10	100	50
262325	65219-238-25	250	30
262450	65219-240-50	500	20

Additional Products & Accessories



freeflex® IV Bags

All the features and advantages of **freeflex+** IV bags, but with a standard injection port.



freeflex® Locking Protective Cap

Bright red cap to indicate drug has been added. Snaps on easily. Offers protection against tampering and discourages diversion.



freeflex® Spiked Secondary Sets

Secure connection for specific use with the **freeflex** IV bags to help avoid spills.



freeflex®+ Vial Adapter with Luer Lock

A convenient twist-on vial adapter for 20mm vials for use with **freeflex+** IV bags.



From **Fresenius Kabi**, one of the largest manufacturers of IV bags in the world.

Fresenius Kabi may be new to you, but our heritage goes back more than 100 years. As a global health care company, we are a leader in high-quality injectable drugs, infusion, transfusion and clinical nutrition products. Headquartered in Europe, we have unique scientific capabilities and state-of-the-art manufacturing and distribution facilities around the globe and across the U.S. We produce more than 1.4 billion units of IV fluids annually.²

Your **ally** in transforming infusion delivery



IVENIX™
INFUSION SYSTEM

Engineered from the ground up to leverage the power of the digital age, Ivenix system is designed with safety in mind.



agilia
Connect

Agilia Connect integrated system of pumps, software, sets, accessories and services.



freeflex®
IV Bags & Solutions

A growing portfolio of IV solutions in **freeflex** and **freeflex+** IV bags.



For more information or to place an order, contact your Fresenius Kabi IV Therapy Sales Representative or call Customer Service at 1.888.386.1300 | www.freeflexivbags.com

2. Fresenius Capital Markets Day 2023 Presentation; Page 25



**FRESENIUS
KABI**

Fresenius Kabi USA, LLC
Three Corporate Drive
Lake Zurich, IL 60047
Phone: 1.888.386.1300
www.fresenius-kabi.com/us

©2024 Fresenius Kabi USA, LLC. All Rights Reserved.
2847-FFX-02-06/22 v3.0