

Drug Reconstitution & Admixture



freeflex®



Trust in Our Experience

Fresenius Kabi is a global health care leader with more than 100 years of experience in injectable drugs, IV solutions, transfusion technologies and nutrition products.

While our IV therapy products have been made in our global network for decades, we recently added a fully automated, vertically integrated plant in Wilson, NC - to help ensure supply reliability, sustainability and responsiveness in the U.S.




Drug Reconstitution and Admixture solutions for the compounding and administration of drugs.

Experience the Difference with **freeflex®+**

From simple needleless drug additions using a Luer lock syringe to whole vials using our single-use **freeflex®+** Transfer Adapters, Fresenius Kabi reconstitution and admixture technologies can help minimize the complexity for pharmacy to compound medications or help make it possible for nursing to admix when needed.

- 1 **freeflex+ IV Bag with Luer Lock Port** - For needleless drug addition up to five times using a Luer lock syringe or using the single-use **freeflex+** Transfer Adapter with 20mm vials.
- 2 **Not Made with PVC or DEHP** - One IV bag platform can be used across a facility for a broad clinical application.
- 3 **Administration Port with Integrated Flow Stop** - Helps prevent free flow from IV bag upon unexpected removal of spike from administration port.
- 4 **Leak Prevention** - Both **freeflex+** ports feature a self-sealing membrane to help prevent leaks after removal of syringe or spike, potentially reducing risk of hazardous drug exposure.
- 5 **Sterility Protection** - Both ports contain easy break-off caps, providing visual tamper evidence and to help ensure sterility prior to first use.
- 6 **Stick Prevention** - Semi-rigid ports may help prevent sharps stick injury during spike insertion.
- 7 **Peel-open Film Overwrap** - Clear, easy-to-open overwrap allows for bag inspection prior to opening.



Swabbable Luer Lock Port



Attaches to a needle-free syringe with one twist

freeflex+, a True Innovation in IV Solutions.

Added flexibility for drug admixture and reconstitution via needleless syringe or **freeflex+** Transfer Adapter.

All the features and advantages of freeflex - plus a unique Luer lock injection port to facilitate needleless drug addition.



For IV fluid management, also check out our standard freeflex IV bags.

For completely needle-free drug addition of up to five 20mm vials, use **freeflex+** Transfer Adapters.



A convenient twist-on **freeflex+** Transfer Adapter for single use with a 20mm drug vial and a **freeflex+** IV bag.



The **freeflex+** portfolio of products

Eight product codes to help meet your reconstitution and admixture clinical needs.

0.9% Sodium Chloride Injection, USP



Product Code	NDC Code	Fill Size (mL)	Expiry	Units per Case
232150	65219-218-50	50	18-Month	60
232200	65219-220-10	100	24-Month	50
232225	65219-222-25	250	36-Month	30
232250	65219-224-50	500	36-Month	20

5% Dextrose Injection, USP



Product Code	NDC Code	Fill Size (mL)	Expiry	Units per Case
252300	65219-234-10	50	18-Month	60
262300	65219-236-10	100	24-Month	50
262325	65219-238-25	250	36-Month	30
262450	65219-240-50	500	36-Month	20

Together, We Can Care Better

Get to know our full line of IV Therapy products.



An expanding line of irrigation solutions featuring **soluflex** irrigation bags and irrigation sets.



Fresenius Kabi offers a full line of IV therapy **hardware, consumables** and **associated accessories**.



Engineered from the ground up to leverage the power of the digital age, **Ivenix** system is designed with safety in mind.



Agilia Connect integrated system of pumps, software, sets, accessories and services.



Our **parenteral nutrition** products are used to help care for critically ill and chronically ill patients in hospitals, long-term care facilities and at home.



We focus on transforming historic breakthroughs in science and health care into **essential medicines**.



For more information or to place an order, contact your IV Therapy Sales Representative or call Customer Service at 1.888.386.1300 | www.freeflexivbags.com



Fresenius Kabi USA, LLC
Three Corporate Drive
Lake Zurich, IL 60047
Phone: 1.888.386.1300
www.fresenius-kabi.com/us

©2025 Fresenius Kabi USA, LLC. All Rights Reserved.
5626-FFX-02-12/24 v2.0